



SPI Supplies 206 Garfield Avenue, West Chester, PA 19380, USA

Use Instructions for Wafer-Mount™ 562

- 1. Cut the Wafer-MountTM sheet to the desired shape and size, and position on a ceramic or glass mounting plate. Use multiple preforms as required to fill in small cavities or gaps caused by warping of the substrate.
- 2. Using a hot plate or oven, heat the mounting plate to 195-210 °F. Make sure to work in a well-ventilated area, and do not overshoot the flow temperature, otherwise, the adhesive will begin to decompose and degrade in strength.
- 3. Using a weight, apply even pressure to the substrate to remove air bubbles and to ensure that the substrate is parallel to the plate.
- 4. Remove the mounting plate from the heat source and allow it to cool slowly to room temperature until the adhesive is hardened. Cool for 20-30 minutes prior to processing.
- 5. Process the substrate as required, then remove parts by re-heating the mounting block to the flow temperature. Use a tool to slide the substrate off the mounting plate.
- 6. For detailed cleaning procedures, refer to the section describing the 562-S Stripper.

Date: 10/15/2015