











Pin Type Mount

	<p>型番: 01506-BA(10), 01506-MB(100) 型番: 01506L-BA(10), 01506L-MB(100) 型番: 01506P-BA(10), 01506-MB(100)</p>	<p>Table diameter: 12.7 mm Pin diameter: 3.2 mm Pin length: 8.0 mm Surface: Lathe, Luster, Polished finish Compatible SEMs: LEO, Cambridge, Philips, ETEC, FEI</p>
	<p>型番: 01506A-BA(10), 01506A-MB(100)</p>	<p>Table diameter: 25.4 mm Pin diameter: 3.2 mm Pin length: 7.5 mm Surface: Lathe finish Compatible SEM: LEO, Leica, Cambridge, Philips, ETEC, FEI</p>
	<p>型番: 01507-BA(10), 01507-MB(100) 型番: 01507L-BA(10), 01507L-MB(100) 型番: 01507P-BA(10)</p>	<p>Table diameter: 12.7 mm Pin diameter: 3.2 mm Pin length: 15 mm Surface: Lathe, Luster, Polished finish Compatible SEM: AMRAY</p>
	<p>型番: 01508-BA(10), 01508-MB(100) 型番: 01508L-BA(10), 01508L-MB(100) 型番: 01508P-BA(10)</p>	<p>Table diameter: 14.3 mm Pin diameter: 6.4 mm Pin length: 16 mm Surface: Lathe, Luster, Polished finish Compatible SEM: AMRAY 900</p>
	<p>型番: 01508A-BA(10), 01508A-MB(100)</p>	<p>Table diameter: 25.4 mm Pin diameter: 6.4 mm Pin length: 15 mm Surface: Lathe finish Compatible SEMs: AMRAY 900</p>
	<p>型番: 01509-BA(10), 01509-MB(100) 型番: 01509L-BA(10), 01509L-MB(100)</p>	<p>Table diameter: 12.7 mm Pin diameter: 3.2 mm Pin length: 15 mm Surface: Lathe, Luster finish Compatible SEMs: AMRAY</p>
	<p>型番: 01509A-BA(10), 01509A-MB(100)</p>	<p>Table diameter: 25.4 mm Pin diameter: 3.2 mm Pin length: 15 mm Surface: Lathe finish Compatible SEMs: AMRAY</p>

	<p>型番: 01517-BA(10), 01517-MB(100) 型番: 01517L-BA(10), 01517L-MB(100)</p>	<p>Table diameter: 12.7 mm Pin diameter: 3.2 mm Pin length: 5 mm Surface: Lathe, Luster finish Compatible SEMs: LEO</p>
	<p>型番: 01519-BA(10), 01519-MB(100) 型番: 01519L-BA(10), 01519L-MB(100)</p>	<p>Table diameter: 12.7 mm Pin diameter: 3.1 mm Pin length: 8 mm Surface: Lathe, Luster finish Compatible SEMs: Zeiss (slot head)</p>
	<p>型番: 156AL-BA(10), 156AL-MB(100) 型番: 156AP-BA(10), 156AP-MB(100)</p>	<p>Table diameter: 25.4 mm Pin diameter: 3.2 mm Pin length: 7.5 mm Surface: Luster, Polished Compatible SEMs: LEO, Leica, Cambridge, Philips, ETEC, FEI</p>